


Respiratory

Nebulizers

Allegiance



AirLife Aerosol Masks

These soft, high-quality AirLife disposable aerosol masks are made from high-grade resins. The feathered and flash free edges and ports provide cool, lightweight comfort and reduce irritation points the wide elastic band with interlocking tabs prevents roll and provides patient comfort all AirLife masks are designed for anatomical compatibility a full line is available in infant, pediatric and adult sizes 50 masks per case

Infant and Pediatric Masks

BGH001261 Pediatric Mask Short Style

BGH001263 Pediatric Mask, Under-the-chin style

BGH001266 Pediatric Dragon Mask

Adult Masks

BGH0012606 Mask Only



AirLife Brand Misty-Neb medication Nebulizers

These disposable Nebulizers are constructed of high-impact plastic to resist breakage and have a built-in anti spill column to minimize loss of costly medication should the unit be tipped over positive threading between the top and bottom of the bottle minimizes leakage of medication during nebulization medication filling is a simple task with the stabilization benefit of built-in legs offered in a broad line of setups AirLife brand Misty-Neb handheld Nebulizers are available with aerosol masks, tee adapters, oxygen tubing, mouthpieces, flextubes and finger controls

Misty-neb Nebulizers, 50/cs

BGH002033 W/adult aerosol mask and 7' (84'') tubing

BGH002043 W/ Pediatric aerosol mask (under-the-chin style)

Misty-neb Nebulizers with Tess Adapter 001508, 50/cs

BGH002037

W/mouthpiece

BGH002034 w/7' (84'') tubing and mouthpiece
BGH002038 w/7' (84'') tubing, mouthpiece, 6' flextube

SUNRISE MEDICAL



DeVilbiss 5650D

The world's most prescribed compressor/Nebulizer, the DeVilbiss Pulmo-Aide, providing effective aerosol treatments for more than thirty years its durability and ease of operation make it the preferred unit for long-term and high-demand users

SUN5650D Pulmo-Aide



DeVilbiss 800

The DeVilbiss Reusable Jet Nebulizers cut treatment time by almost half, are dishwasher safe, can be used with standard compressor, standard mask or used inline with a ventilator circuit 3 models: 800C is green and purple, 800D is blue

and the 800G features and interrupter

SUN Jet Nebulizer



DeVilbiss 4650D-621

DeVilbiss Disposable Nebulizer feature break-resistant plastic that is easy to assemble, clean and provides the patient with 14 days of use capacity of 6.0cc, 50/pk

SUN4650D-621

Disposable Nebulizer, 50/pk

SUN4650D-620

Disposable Nebulizer, 12/pk



DeVilbiss 3650D

The DeVilbiss Pulmo-Aide LT is economical solution for cost-conscious consumers or light duty user comes with a 5-year warranty and includes a disposable Nebulizer

SUN 3650D Pulmo-Aide



DeVilbiss 4650D

The efficient DeVilbiss PulmoMate provides an unbeatable combination of quality and value ideal for traveling it is compact and lightweight easy to use with its extended air outlet and easy access filter comes with a 5-year warranty and includes disposable Nebulizer

SUN 4650D
PulmoMate

MAERSK MEDICAL



Up-Mist Nebulizer

Up-Mist Nebulizer with pediatric aerosol mask and 7 foot oxygen tubing

HOS228 50/cs



UP-Mist Nebulizer

Up-Mist Nebulizer with "T" piece, mouthpiece and 7 foot oxygen tubing

HOS3953 50/cs

Hudson RCI



UP-DRAFT Nebulizer

Self-sealing, high-output Nebulizer available in a variety of configurations produces particles in a therapeutic range continuous flow configuration provides ease of operation, 5-15cc capacity available with 6" reservoir tube all configurations individually packaged, 50/cs

HUD1720 With Tee, Mouthpiece, 7' Tubing

HUD1713 With Pediatric Mask, 7' Tubing



UP-DRAFT II Nebulizer

Smaller sized, high-output medication Nebulizer produces concentrated mist at a flow setting of 6 LPM

Operates even when tilted at 45° angle, 8cc capacity available with 6" reservoir tube anti-spill design prevents loss of medication in any position cap and jar snap apart for easy cleaning available in a variety of individually prepackaged configurations, 50/cs

HUD1732 With Tee, Mouthpiece, 7" tubing
 HUD1734 With Tee, Mouthpiece, reservoir Tubing, 7" Tubing
 HUD1707 With Pediatric Mask, 7" Tubing



MICRO MIST Small Volume Nebulizer

Designed for performance and economy it can be used for hand-held or in-line treatments the MICRO MIST includes the design features needed by today's respiratory care practitioners optimal particle size delivered at efficient nebulization rates with minimal residual consistent performance at angles up to 90° allowing clinicians to administer medication to patients positioned horizontally or with

special treatment requirements easy-seal, threaded cap and 6cc capacity anti-spill jar to prevent leakage available with 6" reservoir tube all configurations individually packaged, 50/cs

HUD1882 With Tee, Mouthpiece, Standard Connector, 7' tubing
 HUD1883 With Tee, Mouthpiece, Reservoir Tubing, standard connector, 7' Tubing
 HUD1886 Pediatric Mask, Standard Connector, 7' Tubing



Pediatric Aerosol Mask

Made from clear soft vinyl for patient comfort adjustable nose clip assures comfortable fit, specifically designed for aerosol therapy connector accepts 22mm aerosol tubing. Tubing not included

HUD1080 Pediatric Aerosol Mask

NELLCOR



Raindrop medication Nebulizer

The Raindrop medication Nebulizer achieves optimum lung penetration for maximum therapeutic benefit documented MMAD particle size range is 1.04-1.10u at 6-8 LPM latex-free components usable with compressor or air outlet, for use with a variety of medications in aerosol therapy single patient use

PUR6-001124-00 w/ Pediatric Aerosol Mask and 7ft Tubing 50/cs
 PUR6-001140-00 w/ Mouthpiece, Tee Adapter, Reservoir and 7ft Tubing, 50/cs

OXYGEN ACCESSORIES

Allegiance



Oxygen Tubing

We offer a very extensive line of disposable vinyl-tipped oxygen tubing it is available with either Smooth-Bor or crush-


resistant lumen to resist occlusion choose from a variety of lengths ranging from 7' (84") to 50' (600")

Tubing with Smooth-Bor lumen

BGH001301 7' (84") 50/cs

Tubing with Crush-Resistant Lumen

BGH001302 7' (84"), 50/cs



AirLife Standard Nasal Cannulas

The AirLife Cannulas offer an angulated, flexible lip plate that provides comfort at the tip of the nasal region, with beveled tips to provide a softer edge. Their thicker tip walls mean good concentricity, furnishing uniform air transport. Disposable Cannulas are curved and are available with flared or nonflared tips, packaged separately or with varying lengths of oxygen tubing.

Cannulas, 50/cs

BGH001310 Flared tip, with 7' (84") Smooth-Bor lumen tubing

BGH001312 Nonflared tip, with 7' (84") Smooth-Bor lumen tubing

Cannulas with Crush-Resistant Lumen Tubing, 50/cs

BGH001321 Flared tip, with 14' (168") tubing

BGH001325 Nonflared tip, with 7' (84") tubing

BGH001326 Nonflared tip, with 14' (168") tubing

AirLife Brand U/Connect-It Female Universal Connectors

The U/Connect-It female universal connector provides a flexible, simple and secure connection to an oxygen Flowmeter or various lengths of AirLife oxygen tubing. There is no need for additional oxygen tubing connectors and Flowmeter adapters. AirLife oxygen tubing, masks, cannulas and Nebulizers features the U/Connect-It female universal connector. Tubing is 7' (84") length, 50 units per case.

U/Connect-It Tubing
BGH001350 Crush-resistant, vinyl-tipped tubing



Tubing Connectors


Lightweight, autoclavable plastic tubing connectors come in a variety of configurations and are transparent to allow easy

observation of flow available sterile in individual packages or clean in dispenser cartons

Clean Connectors

BGH350A 5-in-1 connector, 50/bx, 500/cs

Hudson RCI



Over-the-Ear Cannulas

Very cost and offered in a wide range of configurations, flexible, contoured lip tab provides a high level of patient comfort. Over-the-ear design provides comfortable and secure fit offered with standard (non-flared) or flared nasal tips available with standard or universal connector and standard or crush-proof STAR LUMEN tubing in a variety of lengths, 50/cs.

HUD1100 Adult, Non-Flared Tips, Elastic Head Strap

HUD1101 Pediatric, Non-Flared Tips, Elastic Head Strap

HUD1103 Adult, Standard Tips, 7' Standard Tube

HUD1104 Adult, Flared Tips, 7' Standard Tube

HUD1110 Adult,
Standard Tips, 7' STAR
LUMEN Tubing
HUD1820 Adult,
Softech Design, 7' STAR
LUMEN Tubing

Baxter

Adult Isoflex Breathing Circuit

Parallel Y Bacterial Filter
30" – 90" 20/cs

BGT2L3916

OXYGEN KITS

MADA



Portable Oxygen Cylinder

Portable oxygen cylinder
(685 liters), 1443S
adjustable flow (2-8
LPM) mini-regulator with
extended yoke for use
with humidifier, chrome
cart and nasal cannula
size and color may vary

MAD1630 Port, E-
Cylinder

OXYGEN MASKS

Allegiance



AirLife Oxygen Masks

All AirLife disposable oxygen masks are latex-free and made from soft, high-grade vinyl resins. The feathered and flash-free edges and ports make for cool, lightweight comfort while reducing irritation points. The wide elastic band with interlocking tabs and no metal clips prevents roll and provides patient comfort. All AirLife oxygen masks are designed for anatomical compatibility. A full range is available to meet all requirements in infant, pediatric and adult sizes. The AirLife 3-in-1 masks 001202, 001203, 001362 and 001210 are actually three masks in one. Masks are packaged as high concentration. Nonbreathing masks complete with reservoir bags and one-way valves, 50 masks/cs. **Pediatric Vinyl Oxygen Masks (under-the-chin-style)**

BGH001262 Medium
Concentration mask with
7' (84") tubing

Pediatric Vinyl Oxygen Masks (Short)

BGH001260 Medium
Concentration Mask with
7' (84") tubing

Adult Vinyl Oxygen Masks (under-the-chin- style)

BGH001201 Medium
concentration mask with
7' (84") tubing

BGH001203 High
concentration,
nonbreather 3-in-1
mask with safety vent and
7'


Hudson RCI



Oxygen Masks

Made from clear, soft
vinyl for patient comfort
and visual patient
assessment. Adjustable
nose clip assures
comfortable fit. Complete
with 7-ft oxygen supply
tubing. Individually
packed, 50/cs.

HUD1035 Child,
Medium Concentration
HUD1040 Adult,
Medium Concentration
HUD1041 Adult,
Elongated, Medium
Concentration



Non-Rebreathing Mask

Made from clear, soft vinyl for patient comfort and visual patient assessment adjustable nose clip assures comfortable fit low-resistance check valve prevents Rebreathing and allows exhaled gas to escape complete with 7-ft oxygen supply tubing individually packaged, 50/cs

HUD1060 Adult

Non-Rebreathing mask with Safety Vent

Low-resistance check valve prevents Rebreathing and allows exhaled gas to escape safety vent allows for entrainment of room air complete with 7-ft oxygen supply tubing adjustable nose clip assures comfortable fit individually packaged, 50/cs

HUD1059 Adult

PULMONARY FUNCTION

Hudson RCI



POCKETPEAK Peak Flow Meter

Compact and easy to use, and now redesigned with durable, sliding zone indicators which encourage patient compliance colorful markers allow patients to easily monitor changes in airflow, and healthcare professionals to recommend a treatment plan of action when readings fall in the red, yellow or green zones meets or exceeds all National Asthma Education Program and National Heart, Lung, Blood Institute standards ABS plastic housing and tempered stainless steel vane are completely dishwasher-safe the reusable mouthpiece attaches to the POCKETPEAK for secure storage Scale increments: 10 LPM

HUD1801
POCKETPEAK

POCKETPEAK Mouthpiece

Disposable cardboard mouthpiece packaged in convenient dispenser box, 100/bx

HUD1805 Disposable cardboard mouthpiece

MADA



Assess Peak Flow Meter

Standard adult range (60-880 LPM)

MAD710 Peak Flow Meter, 12/cs

MAD711 Disposable Mouthpiece, 100/bx

Puritan Bennett



PocketPeak Peak Flow Meter

Helps patients monitor their asthma easily and affordably and promotes compliance to daily peak flow monitoring one-size-fits-all shape is ergonomically right for patients of all ages and lungs disorders scales suit

most age groups not positional dependent meets standards of the American Thoracic Society and the National Asthma Education and prevention program for accuracy and reliability stick-on zone labels for individualized asthma management Readings are simple to interpret and understand small and lightweight to fit pocket or purse instruction manual/patient diary is easy to read and understand mouthpiece fits adult or pediatric slidable zone scales for easy adjustment
PURP-063436-00
PocketPeak

Flow Sensors
PURP-000263-00
FSII Flow Sensors,
50/pkg

Spirometrics



Spir-O-Flow Peak Flow Pocket Monitor

The Spir-O-Flow features a dual scale that allows for both standard and low range measurements, is pocket size, dishwasher safe, and meets the National Asthma Education Program technical standard

packaging comes with instructions, mouthpiece, and daily record chart in a zip-lock pouch

SPI2350 Spir-O-Flow Peak Flow Pocket Monitor

Accessories

SPI2175 Pulmonary Filter and Mouthpiece
SPI2014 Adult Mouthpiece, Plastic
SPI3014 Pediatric Mouthpiece, Cardboard
SPI2104 Noseclip, 10/bx
SPI143 Filter Disk, 50/pkg
SP2128 printer Paper, Flowmate

Trach Trays & Supplies

Hudson RCI



SOFTECH

Endotracheal Tubes
SOFTECH oral-nasal endotracheal tubes with Murphy eye combine quality, dependability and performance high volume/low pressure dipped cuff provides

secure, a traumatic seal ergonomically designed pilot balloon with large surface area provides maximum tactile feedback patented beveled 15mm connector allows easy reinsertion after tube is trimmed tubes individually packaged, sterile, 10/bx

HUD340070 7mm

SHERIDAN/HVT Tracheal Tubes

High volume, large diameter cuff provides positive tracheal wall seal and reduces area of tracheal wall contact low pressure cuff minimizes capillary restriction for oral or nasal use available with Murphy eye only tubes individually packaged, sterile, 10/bk

RUSH



Endotracheal Tubes

Cuffed Murphy Eye
Sterile, oral or nasal use

RUS112082050
5.0mm, 10/bx
RUS112082065
6.5mm, 10/bx
RUS112082070
7.0mm, 10/bx
RUS112082075
7.5mm, 10/bx

RUS112082080
8.0mm, 10/bx

DALE



DALE Tracheostomy Tube Holder

Latex-Free to make trach tube stabilization more secure and to avoid the knotty problems associated with the tying and untying of twill ties, consider the Dale Tracheostomy Tube Holder the narrow fastener tabs of the Dale Tracheostomy Tube Holder engage and disengage the trach plate quickly and easily, so tube stabilization is faster and safer than ever before and unlike unyielding twill tape, the soft neckband of the Dale Tracheostomy Tube Holder incorporates a length of stretch material that allows for cough reflex the stretch material also accommodates any edema around the neck, and ensures a safe, snug and comfortable fit an extension piece is available to accommodate your bariatric patients

DAL240 Trach
Tube Holder Strap